



DISPOSABLE BIFLEX GOWN



Product Description

Developed in Opia laboratories, our Biflex Surgical Gowns are made with spunbond nonwoven interlining, which is highly absorbent with hydrophilic application; It is produced for disposable operating room drape-compress/drape laminated with hot melt (hot melt) of microemboss engraved barrier film. It is produced in 2 layers. Thanks to its high absorbency, our medical fabrics provide protection by absorbing the liquid very quickly and keeping it in the barrier layer.

Our fabrics produced according to EN 13795 standards; It is produced at least 2 times above the minimum limit values to be reached in the operating room environment; It performs above the EN 13795 test criteria.

BIFLEX GOWN STANDARD

It is breathable, waterproof, lightweight and easy to wear. Knitted or elastic cuffs are used optionally. Our products can be easily used by patients' relatives, patients, nurses, doctors and laboratory workers. It is produced as a double-tie waist and neck-tied or velcro system. Our products are produced as Sterile and Non-Sterile.

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

ORDER REF. #	OP-00205
WEIGHT	40-50 gsm
COLOURS AVAILABLE:	Medical Blue, Blue, Navy Blue, Green, White
SIZES AVAILABLE:	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	HOSPITAL, LABORATUVARY
MATERIAL:	PP+SPUNLACE (NONWOVEN)
PER CASE:	60 STERILE - 80 (NONSTERIL) PIECES
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET BOX:	24



MDR **CE** 2696

CERTIFICATION NO: M.2022.MDR.1002

THIS PRODUCT:

UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA
NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO
THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN
PARLIAMENT. COMPLIED WITH THESE STANDARDS:
CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,
EKOTEKS: AAMI PB-70, AATCC -42, AATCC-127, LEVEL II-III-IV



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY