



INDUSTRIAL TYPE 5/6 OVERALLS



Product Description

As Opia Medikal & Health, MDR Type 5 and Type 6 disposable overalls, produced in Protective PPE, MDR standards, provide maximum protection without sacrificing comfort. It also has a shelf life of 5 years. Our industrial Type 5 and Type 6 protective overalls are manufactured in our factory specially designed in accordance with European standards. Disposable coveralls provide a high level of protection for your comfort with both Microporous (High Protection) and SMMS material, as well as a special lamination that provides maximum protection while keeping you cool. Low cost hazardous dust protection The Opia™ Protective Coverall is an economical solution for protection from hazardous dusts such as asbestos fibers exposed during site inspection. While protecting the outside, it absorbs the sweat of the employee thanks to its specially produced fabric. Besides, the breathable material allows air to escape inside to reduce heat build-up and discomfort. Available in white, blue, red and orange, the series is CE Type 5 and Type 6 approved for both hazardous dusts and liquid splashes.

Scope of application

Asbestos inspection • Site preparation, insulation • Construction • Cement manufacturing • Grinding and polishing • Light industrial cleaning • Machine maintenance. Our products have MDR Certificate, FDA Certificate and GMP Certificates. In addition, it is produced according to the national medical production license.

Technical Standards

Material : SMMS + Polypropylene + Microporous film coated

Appearance: Hooded, elasticated wrists, waist and ankles.

Strength : Excellent resistance to liquid and chemical splash

Impermeability: All parts that are stitched are closed with a strip by applying a special application at certain temperatures.

Front Closure: Zipper system and on it, zipper closure according to the size of each product

Our products are produced as Sterile and Non Sterile.



EN 14126 EN 1073-2 EN 1149-5 EN 13982-1 EN 13034 EN 340 EN 14605

MDR CE 2696 **CERTIFIED PRODUCT GUARANTEE**

CERTIFICATION NO : M.2022.MDR.1002

THIS PRODUCT:
 UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA NOTIFICATION BOARD:
 2696 CERTIFIED IN APRIL 2022 ACCORDING TO THE 2017/745 REGULATION OF THE COUNCIL OF
 THE EUROPEAN PARLIAMENT, COMPLIED WITH THESE STANDARDS:
 CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,
 EKOTEKS : AAMI PB-70 ,AATCC -42 ,AATCC-127,LEVEL II-III-IV

SIZE	LENGTH	WIDTH	ARM LENGTH
S	185	155	58
M	190	160	60
L	193	163	63
XL	200	167	66
2XL	203	170	68
3XL	207	173	70

ORDER REF. #	END TYPE 5-603 / END TYPE 6-604
WEIGHT	55-60 GSM
COLOURS	White, Red, Blue, Yellow
SIZES	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	In Hospitals, Factories, Chemical Manufacturing
MATERIAL:	SMMS + Polypropylene + Microporous film coated
IN BOX	50 PIECES
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET	24 KUTU



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
 No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY