



# MDR TYPE 3/4 OVERALLS



## Product Description

Opia Medikal & Health, our disposable MDR - Type 3 and Type 4 coveralls are highly resistant to microporous, dry particles and liquid chemical splashes. Our laminated microporous Nonwoven fabrics and MDR Type 3 and Type 4 coveralls have a breathable structure. It is extremely comfortable for long-term use.

All the necessary standards are applied over the FDA / PPE Regulations in the processes until the final stage of our MDR TYPE 3 and TYPE 4 overalls, which starts with the technical testing of the raw materials before production. The most important part of our products, which go through many stages in production, is keeping the sealing at the top level. The combination of microporous Coverall, soft polypropylene non-woven fabric and microporous film allows moisture vapor to escape to keep the wearer comfortable. It has good barrier property for wet or liquid and dry particles. It provides good protection in highly sensitive environments including medical applications, pharmaceutical factories, sterile environments, non-toxic liquid handling operations and general industrial work areas. Ideal for Security, Mining, Clean room, Food industry, Medical, Laboratory, Pharmaceutical, Industrial pest control, Machine Maintenance and Agriculture. In addition, each of our products that are produced are checked one by one by the quality control unit and packaged with the "quality controlled" stamp. Our products have MDR Certificate and GMP Certificates. In addition, it is produced according to the national medical production license.

## Technical Standards

Material : SMS + Polypropylene + Microporous film coated

Appearance: Hooded, elasticated wrists, waist and ankles.

Strength : Excellent resistance to liquid and chemical splash

Impermeability: All parts that are stitched are closed with Blue Stripe by applying specially at certain temperatures.

Front Closure: Zipper system and on it, zipper closure according to the size of each product

SIZE	LENGTH	WIDTH	ARM LENGTH
S	185	155	58
M	190	160	60
L	193	163	63
XL	200	167	66
2XL	203	170	68
3XL	207	173	70

ORDER REF. #	MDR TYPE 3-501 / MDR TYPE 4-502
WEIGHT	55-60 GSM
COLOURS	White, Red, Blue, Yellow
SIZES	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	In Hospitals, Factories, Chemical Manufacturing
MATERIAL:	SMMS + Polypropylene + Microporous film coated
IN BOX	50 PIECES
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET	24 KUTU



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepke Plaza  
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY