



DISPOSABLE CHEMOTHERAPY GOWN



Product Description

Our OPIA Chemotherapy gowns are an advanced gown with excellent protection against acids, splashes and high-risk contaminants. It is most commonly used in chemical laboratories, high-risk drug or drug manufacturing facilities, hospital environments, and hazardous manufacturing facilities. It provides user comfort and protection with its latex-free and nitrile-free tested fabric. OPIA Chemotherapy Gowns are meticulously crafted from a soft, multi-layered non-woven spunbond material with a liquid-proof polyethylene coating for exceptional liquid and chemical resistance. It has ultrasonically sealed seams for superior protection and high tensile strength.

These aprons also feature a highly visible blue band to ensure the right apron is worn in the right environment. Velcro closure at the neck wraps the apron, making it easy to remove the set. Extra long laces combined with multiple size combinations provide a comfortable and secure fit. Extra long ties also tie the knot to the wearer's side, allowing the suit to be secured outside of compound activities.

CHEMOTHERAPY GOWN STANDARD

Disposable, lint-free and resistant to chemotherapy permeability long sleeves and knit cuffs, hook-and-loop fastening at the back front closed

It is meticulously produced from polyethylene-coated and multi-layered nonwoven spun bond material. Meets USP <800> standards for handling hazardous drugs.

To ensure the safe use of OPIA Chemotherapy gowns, personnel should wear them only during patient care and interaction, or during activities that require protection such as handling, preparation, post-treatment care, and disposal to ensure safe use. Personnel should also check the apron for visible defects such as defects or tears before putting it on and replace it after spills, splashes or 2-3 hours of continuous use.

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

ORDER REF. #	OP-00206
WEIGHT	45 gsm + 15 PE
COLOURS AVAILABLE:	Medical Blue, Blue, Navy Blue, Green, White
SIZES AVAILABLE:	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	In chemistry laboratories, high-risk drug or drug manufacturing facilities, hospital environments and hazardous production facilities.
MATERIAL:	SMS +PE
PER CASE:	Sterile 45 Pieces /Nonsterile 60 Pieces
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET BOX:	24



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY