



PP+PE GOWN



Product Description

OPIA PP+PE Disposable Surgical Gown is made of high quality 25gr/m2 PP+PE fabric and it offers great comfort while using the products. it is suitable for daily use in hospitals and it can be easily worn over a daily dress. Optionally, the products are available with or without strips. Products it can be used by doctors and patients. it has 2 waist ties. Neck lacing and velcro closure. Our PP+PE Disposable surgical gowns come in EO sterile packaging with the product label indicating the products. it can also be non-sterile on request.

- * PP+PE Fabric (Polypropylene Nonwoven Fabric + Polyethylene Nonwoven Fabric)
 - * Non Sterile (Autoclavable at 1340)
 - * Anti-Alcohol, Anti-Blood, Anti-Oil, Anti-bacteria
 - * Waterproof, Alkali Resistant
 - * Produced from 25 gr/m2 PP+PE Fabric.
 - * Sleeves are available in both strip and rubber versions it can be easily worn over a dress.
 - * Breathable, Flexible
 - * Intended Use: it is for the protection of the body against drop particles and infected agents coming from outside.
- Advantage:
- * Durable
 - * Available in striped or unstriped design *Easy to wear
 - * Breathable and Comfortable *High quality insulation

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

ORDER REF. #	OP-00207
WEIGHT	40 gsm-55 gsm
COLOURS AVAILABLE:	Medical Blue, Blue, White
SIZES AVAILABLE:	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	Hospitals, Chemical Factories, Petroleum Refineries, Construction Industry Etc.
MATERIAL:	PP+PE (Polipropilen + Polietilen)
PER CASE:	60 PIECES
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET BOX:	24








MDR CE 2696

CERTIFICATION NO:M.2022.MDR.1002

THIS PRODUCT:
 UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA
 NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO
 THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN
 PARLIAMENT. COMPLIED WITH THESE STANDARDS:
 CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,
 EKOTEKS : AAMI PB-70, AATCC -42, AATCC-127, LEVEL II-III-IV






PHONE : +90 (850) 308 6742
 WHATSAPP : +90 (533) 636 6742
 Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
 No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY