



# DISPOSABLE UROLOGY GOWN



## Product Description

OPIA UROLOGY Apron is for single use only. Hybrid Fabric PP Spunbond+PE Film(Biflex) Ultrasonically welded, ergonomically designed to provide a comfortable sitting position for surgeons. it has excellent liquid impermeability properties to protect the operating room personnel from possible contamination. Thanks to its nonwoven fabric Prevents cross contamination by reducing pilling. it is resistant to liquids with its barrier feature. it is extremely resistant to moisture and bacterial penetration. Our UROLOGY gowns, which are resistant to puncture, tear and abrasion, do not contain toxic substances. Apron cuff is compatible with surgical gloves, prevents gloves from slipping. Our fabrics produced according to EN 13795 standards; it is produced at least 2 times above the minimum limit values to be reached in the operating room environment; Performs above EN 13795 test criteria.

### EO Sterilize

Hybrid Fabric (SMS & PP Spunbond + PE Film (BIFLEX)  
Ultrasonic Welded

Designed to provide comfortable sitting position to surgeon  
Surgically Folded

Excellent water repellence to protect surgeons from contamination. Prevention of cross contamination by reducing lint  
Resistant to liquids with its barrier feature. Resistant to moisture and bacterial penetration. Puncture, tear and abrasion resistant.

Free from toxic substances

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

<b>ORDER REF. #</b>	OP-00204
<b>WEIGHT</b>	60 gsm
<b>COLOURS AVAILABLE:</b>	Medical Blue, Green
<b>SIZES AVAILABLE:</b>	S, M, L, XL, 2XL, 3XL
<b>TYPICAL APPLICATIONS:</b>	HOSPITAL
<b>MATERIAL:</b>	SMS + PP SPUNBOND + PE (BIFLEX)
<b>PER CASE:</b>	Sterile 50 Pieces /Nonsterile 80 Pieces
<b>BOX SIZE:</b>	40X60X40
<b>PALLET SIZE:</b>	80X120
<b>PALLET BOX:</b>	24



**MDR** **2696**  
**CERTIFICATION NO:M.2022.MDR.1002**

THIS PRODUCT:  
UDM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA  
NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO  
THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN  
PARLIAMENT. COMPLIED WITH THESE STANDARDS:  
CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,  
EKOTEKS: AAMI PB-70, AATCC -42, AATCC-127, LEVEL II-III-IV



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza  
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY