



# DISPOSABLE STANDARD GOWN



## Product Description

OPIA disposable standard surgical gowns are designed to reduce risk and adapt to short-term operations. It ensures that the healthcare professional, who does not have direct contact with the surgical field in high-risk operations, remains sterile and works comfortably. OPIA disposable standard surgical gowns, produced from OPIA SMS fabrics, are designed to keep the user's protection and comfort at the highest level. The long and elastic cuffs of disposable surgical gowns also eliminate the risk of the glove slipping during the operation. OPIA disposable standard gowns are lightweight SMS gowns made with an additional liquid repellent treatment to increase protection. SMMS technology provides high comfort and breathability for an effective gown. If OPIA disposable standard surgical gowns are sterile, they are packed in sterilization bags. OPIA disposable standard surgical gowns can be packed with medical towels depending on the user's preference. It can also be wrapped in a pack sleeve. Our disposable standard surgical gowns are available in a wide variety of sizes.

Disposable surgical gown

individually sterile packaged

Made of liquid-repellent fleece with a micro-fiber layer

Wear resistant and low lint

High resistance to moisture penetration and migration of microbes, Fixed fastening straps and ultrasonic or overlock seams, Elastic cuffs on sleeves Individually adjustable hook and loop closure neck

OPIA Standard surgical gowns are available in sizes S-3XL.

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

<b>ORDER REF. #</b>	OP-00201
<b>WEIGHT</b>	35/40/45/50 gsm
<b>COLOURS AVAILABLE:</b>	Medical Blue, Blue
<b>SIZES AVAILABLE:</b>	S, M, L, XL, 2XL, 3XL
<b>TYPICAL APPLICATIONS:</b>	In Hospitals, Healthcare Facilities
<b>MATERIAL:</b>	SMS (Spunbond+Meltblown+Spunbond)
<b>PER CASE:</b>	Sterile 60 Pieces /Nonsterile 80 Pieces
<b>BOX SIZE:</b>	40X60X40
<b>PALLET SIZE:</b>	80X120
<b>PALLET BOX:</b>	24



**MDR CE 2696**

**CERTIFICATION NO:M.2022.MDR.1002**

THIS PRODUCT:  
UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA  
NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO  
THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN  
PARLIAMENT. COMPLIED WITH THESE STANDARDS:  
CENTROCOT ISO 13795-1, ISO 22612-4 ,AATCC 9372  
EKOTEKS : AAMI PB-70 ,AATCC -42 ,AATCC-127,LEVEL II-III-IV



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza  
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY