



# DISPOSABLE ISOLATION GOWN



## Product Description

Our Opia Medical Disposable Isolation Gown is designed to prevent infectious diseases and pathogens from infecting your clothes and spreading them to others. There are two types of isolation gowns; liquid resistant isolation gowns and impermeable isolation gowns. Liquid-resistant gowns repel light liquids such as blood and urine, but only for a short time as these can be absorbed into the gown after a while. Impermeable isolation gowns can withstand greater amounts of liquid and will not seep through the material and contaminate the garment. It is worn to protect healthcare workers, use for visitors, when the risk of fluid contact or penetration is low, including blood collection, suturing, or while working in a pathology laboratory. Production is made in various colors and sizes.

### Technical Standards

It is breathable, waterproof, lightweight and easy to wear.

Knitted or elastic cuffs are used optionally.

Our products can be easily used by patients' relatives, patients, nurses, doctors and laboratory workers.

It is produced as a double-tie waist and neck-tied or velcro system.

Our products are produced as Sterile and Non-Sterile.

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

<b>ORDER REF. #</b>	OP-00101
<b>WEIGHT</b>	20/25/30/40 gsm
<b>COLOURS AVAILABLE:</b>	Medical Blue, Blue, Navy Blue, Yellow, White
<b>SIZES AVAILABLE:</b>	S, M, L, XL, 2XL, 3XL
<b>TYPICAL APPLICATIONS:</b>	Laboratory, dust-free, workshop, medical care, electronic factory, food industry, beauty care, decoration engineering, hardware industry, cleaning and daily life, etc.
<b>MATERIAL:</b>	SS (Spunbond + Spunbond)
<b>PER CASE:</b>	120 Piece
<b>BOX SIZE:</b>	40X60X40
<b>PALLET SIZE:</b>	80X120
<b>PALLET BOX:</b>	24



**MDR** **CE** 2696

**CERTIFICATION NO: M.2022.MDR.1002**

THIS PRODUCT:  
UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA  
NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO  
THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN  
PARLIAMENT COMPLIED WITH THESE STANDARDS:  
CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372  
EKOTEKS: AAMI PB-70, AATCC -42, AATCC-127, LEVEL II-III-IV



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza  
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY