



DISPOSABLE REINFORCED GOWN



Product Description

Our OPIA Sterile Reinforced Surgical Gowns are made to provide the operator with high thermal comfort, high flexibility. It is made of hydrophobic soft touch nonwoven material with a minimum weight of 46 gsm SMS and a hydrophilic PE laminated spunbond backsheet in critical areas. Our Sterile Reinforced Surgical Gowns are medically folded to enable aseptic administration. Our Sterile Reinforced Surgical Gown is designed to be used to protect the healthcare professional during the surgical procedure and to provide the user with absolute comfort in use. It is mostly designed for surgeons. With its hydrophobic feature and hydrophilic reinforcement, the sterile reinforced surgical gown is designed to withstand any liquid from side to side and facilitate the surgeon's movement during surgery. This type of surgical gown is mostly used in standard and long-term surgical operations. EO Sterilized - Extra Reinforced layer in the arms and chest Level III-Level III AAMI - Nonwoven Materials (SMS) Ultrasonic Welded - Surgically Folded Excellent water repellence to protect surgeons from contamination Prevention of cross contamination by reducing lint Resistant to moisture and bacterial penetration Puncture, tear and abrasion resistant. Latex free - Does not contain toxic substances

SIZE	LENGTH	WIDTH	ARM LENGTH
S	110	135	58
M	115	140	59
L	120	145	60
XL	125	150	61
2XL	130	155	62
3XL	135	160	63

ORDER REF. #	OP-00202
WEIGHT	46-51 gsm
COLOURS AVAILABLE:	Medical Blue, Blue
SIZES AVAILABLE:	S, M, L, XL, 2XL, 3XL
TYPICAL APPLICATIONS:	In Hospitals, Healthcare Facilities
MATERIAL:	SMS (Spunbond+Meltblown+Spunbond) + PE
PER CASE:	Sterile 50 Pieces /Nonsterile 65-70 Pieces
BOX SIZE:	40X60X40
PALLET SIZE:	80X120
PALLET BOX:	24



MDR CE 2696
CERTIFICATION NO:M.2022.MDR.1002



THIS PRODUCT:
UDEM Adriatic D.O.O RADNICKA CESTA 54/R3 10000 ZAGREB CROATIA
NOTIFICATION BOARD: 2696 CERTIFIED IN APRIL 2022 ACCORDING TO
THE 2017/745 REGULATION OF THE COUNCIL OF THE EUROPEAN
PARLIAMENT. COMPLIED WITH THESE STANDARDS:
CENTROCOT ISO 13795-1, ISO 22612-4, AATCC 9372,
EKOTEKS : AAMI PB-70, AATCC -42, AATCC-127, LEVEL II-III-IV



PHONE : +90 (850) 308 6742

WHATSAPP : +90 (533) 636 6742

Address : İncilipınar Mah. 36006 Cad.Kepkep Plaza
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY