



TYPE 5/6 OVERALLS



Product Description

As Opia Medikal & Health, TYPE 5 and TYPE 6 disposable overalls, produced in protective PPE, MDR and FDA standards, provide maximum protection without sacrificing comfort. Type 5 and Type 6 protective disposable coveralls provide a high level of protection for your comfort with both Microporous (High Protection) and SMS material and also with a special lamination that provides maximum protection while keeping you cool. As one of the leading suppliers in Turkey, we continue to work with our R&D and Innovation team to ensure that our disposable medical coveralls are durable and recyclable, unlike the cheap and bad disposable coveralls that we see torn when you wear them. In addition, every product that is produced is checked one by one by the quality control unit and packaged with the "quality controlled" stamp. Our products have MDR Certificate, FDA Certificate and GMP Certificates. In addition, it is produced according to the national medical production license.

Technical Standards

Material : Polypropylene + Microporous film coated
Appearance : Hooded, elasticated wrists, waist and ankles.
Strength : Excellent resistance to liquid and chemical splashes
Sealing : All parts that are stitched are closed with stripe by applying special applications at certain temperatures.
Front Closure : Zipper system and on it, zipper closure according to the size of each product
Our products are produced as Sterile and Non-Sterile.

| SIZE | LENGTH | WIDTH | ARM LENGTH |
|------|--------|-------|------------|
| S | 185 | 155 | 58 |
| M | 190 | 160 | 60 |
| L | 193 | 163 | 63 |
| XL | 200 | 167 | 66 |
| 2XL | 203 | 170 | 68 |
| 3XL | 207 | 173 | 70 |

| | |
|------------------------------|---|
| ORDER REF. # | OIM-0301 |
| WEIGHT | 55-60 GSM |
| COLOURS | White, Red, Blue, Yellow |
| SIZES | S, M, L, XL, 2XL, 3XL |
| TYPICAL APPLICATIONS: | In Hospitals, Factories, Chemical Manufacturing |
| MATERIAL: | SMMS + Polypropylene + Microporous film coated |
| IN BOX | 50 PIECES |
| BOX SIZE: | 40X60X40 |
| PALLET SIZE: | 80X120 |
| PALLET | 24 KUTU |



PHONE : +90 (850) 308 6742
WHATSAPP : +90 (533) 636 6742
Address : İncilipınar Mah. 36006 Cad.Kepke Plaza
No: 20/29 Kat: 4 Şehitkamil Gaziantep / TURKEY